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WINT

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- Akiyama H, Oono T, Saito M, Iwatsuki K (2004) Assessment of cadexomer iodine against *Staphylococcus aureus* biofilm in vivo and in vitro using confocal laser scanning microscopy. *J Dermatol* 31(7): 529–34
- Bain K, Hansen AS (2020) Strengthening implementation success using large-scale consensus decision-making - A new approach to creating medical practice guidelines. *Eval Program Plann* 79: 101730
- Dowsett C, Muentert C, Bain K, Bain M (2020) Closing the gap between the evidence and the bedside: A consensus report on Exudate Management. *Wounds International* 11(3): 64–8
- Dowsett C, Swanson T, Karlsmark T (2019) A Focus on the Triangle of Wound Assessment – addressing the gap challenge and identifying suspected biofilm in clinical practice. *Wounds International* 10(3): 16–21
- European Wound Management Association (2005) *Position Document: Identifying Criteria for Wound Infection*. MEP Ltd, London. Available at: http://www.cslr.cz/download/English_pos_doc_final.pdf (accessed 24.11.2020)
- Guest JF, Ayoub N, McIlwraith T et al (2015) Health economic burden that wounds impose on the National Health Service in the UK. *BMJ Open* 5(12): e009283
- International Wound Infection Institute (IWII) Wound infection in clinical practice. *Wounds International* 2016
- Johani K, Malone M, Jensen S et al (2017) Microscopy visualisation confirms multi-species biofilms are ubiquitous in diabetic foot ulcers. *Int Wound J* 14(6): 1160–9
- Keast D, Bain K, Hoffmann C et al (2020) Managing the Gap to promote Healing in Chronic Wounds – An International Consensus. *Wounds International* 11(3): 58–63
- Keast D, Swanson T, Carville K et al (2014) Ten Top Tips: Understanding and managing wound biofilm. *Wounds International* 5(2): 20–4
- Malone M, Bjarnsholt T, McBain AJ et al (2017) The prevalence of biofilms in chronic wounds: a systematic review and meta-analysis of published data. *J Wound Care* 26(1): 20–5
- Moore Z, Strapp H (2015) Managing the problem of excess exudate. *Br J Nurs* 24(15): S12–7
- Moura MRT, Dowsett C, Bain K, Bain M (2020) Advancing practice in holistic wound management: a consensus-based call to action. *Wounds International* 11(4): 70–5
- Murphy C, Atkin L, Swanson T et al (2020) International consensus document. Defying hard-to-heal wounds with an early antibiofilm intervention strategy: wound hygiene. *J Wound Care* 29(Suppl 3b): S1–28
- Ousey K, Atkin L (2013) Optimising the patient journey Made Easy. *Wounds UK* 9(2): 1–6
- Ousey K, Cook L (2012) Wound assessment made easy. *Wounds UK* 8(2): 1–4
- Penhallow K (2005) A review of studies that examine the impact of infection on the normal wound-healing process. *J Wound Care* 14(3): 123–6
- Persival SL (2017) Importance of biofilm formation in surgical infection. *BLS Society* 104(2): 85–94
- Percival SL, Bowler P, Woods EJ (2008) Assessing the effect of an antimicrobial wound dressing on biofilms. *Wound Repair Regen* 16(1): 52–7

- Phillips PL, Yang Q, Davis S, et al (2015) Antimicrobial dressing efficacy against mature *Pseudomonas aeruginosa* biofilm on porcine skin explants. *Int Wound J* 12(4): 469–83
- Phillips PL, Yang Q, Schultz GS (2013) The effect of negative pressure wound therapy with periodic instillation using antimicrobial solutions on *Pseudomonas aeruginosa* biofilm on porcine skin explants. *Int Wound J* 10(1): 48–55
- Phillips PL, Wolcott RD, Fletcher J, Schultz GS (2010) Biofilms Made Easy. *Wounds International* 1(3): 1–6
- Rodeheaver GT, Ratliff CR (2007) Wound cleansing, wound irrigation, wound disinfection. In: Rodeheaver GT et al (eds.) *Chronic Wound Care: A Clinical Source Book for Healthcare Professionals*. HMP Communications, Malvern, PA: pp331–42
- Scali C, Kunimoto B (2013) An update on chronic wounds and the role of biofilms. *J Cutan Med Surg* 17(6): 371–6
- Schultz G, Bjarnsholt T, James GA et al (2017) Consensus guidelines for the identification and treatment of biofilms in chronic nonhealing wounds. *Wound Repair Regen* 25(5): 744–57
- Swanson T, Keast D, Cooper R et al (2015) Ten Top Tips: identification of wound infection in a chronic wound. *Wounds International* 6(2): 22–7
- Swanson T, Grothier L, Schultz G (2014) *Wound Infection Made Easy*. Wounds International, London.
- Wolcott RD, Kendra P, Rumbaugh GJ et al (2010) Biofilm maturity studies indicate sharp debridement opens a time-dependent therapeutic window. *J Wound Care* 19(8): 320–8
- World Union of Wound Healing Societies (2019) *WUWHS Consensus Document: Wound Exudate: Effective Assessment and Management*. Wounds International, London. Available at: <https://bit.ly/3692LbN> (accessed 24.11.2020)
- World Union of Wound Healing Societies (2016) *Position Document: Management of Biofilm*. Wounds International, London. Available at: <https://bit.ly/3l3k1mU> (accessed 24.11.2020)
- World Union of Wound Healing Societies (2008) *Wound Infection in Clinical Practice. A WUWHS International Consensus*. London, MEP Ltd. Available at: <https://bit.ly/39ixa9x>(accessed 24.11.2020)
- Wounds UK (2018) *Best Practice Statement: Improving Holistic Assessment of Chronic Wounds*. Wounds UK, London. Available at: <https://bit.ly/3fBggDY> (accessed 24.11.2020)
- Wounds UK (2017) *Best Practice Statement: Making Day-to-Day Management of Biofilm Simple*. Wounds UK, London. Available at: <https://bit.ly/2JdBnjG> (accessed 24.11.2020)